How to Download NMR data

20200707 lzhu

 First of all, you must connect <u>UCR VPN</u> except your PC is hard wired to Ethernet cable inside Chemical Science Building. Please see this link about how to use UCR VPN: <u>https://ucrsupport.service-</u>

now.com/ucr portal/?id=kb article&sys id=0bf268b81bfd9490f3444158dc4bcbb3

2. You will need to use SFTP tool to download data, e.g. WinSCP.

All the NMR data servers have their own IP addresses: AV700: acif700.dyn.ucr.edu AV600: acif600.dyn.ucr.edu NEO600: acifneo600.dyn.ucr.edu AV500: acif500.dyn.ucr.edu NEO400: acifneo400.dyn.ucr.edu

AV300: acif300.dyn.ucr.edu

Username (usually your PI's last name) and Password are the same when you use the corresponding NMR PC. Use the **Port: 22** when you connect.

The following photo is an example connection in WinSCP:

nmrsu - av600 - WinSCP		- 🗆 X
Local Mark Files Commands Session Optio	ns Remote Help	
🕀 📚 Şynchronize 🗾 🦑 💽 🔅	😭 Queue 🗸	Transfer Settings Default 🔹 🖉 🗸
av600 × 🗳 New Session		
🛄 Desktop 🔹 🚰 🕶 🔽 🖛 🔹 🛶 🔹 🙀	a 🗈 🏠 🎜 🖏	🔄 nr 🔹 🚰 🔹 🟹 🔹 🖛 🔹 🖚 🖈 💼 🔂 🏠 🎜
	en lienen beschlichense prinder	📲 Download 🗸 📝 Edit 🤟 🗱 🔂 🕞 Properties 🛛 🔭 🛛 🛨
C:\Users\Danb\Desktop\		/home/opt/nmrsu/
Name User List UCR Your PC disk Temp Software Icons Personal NMR STD Bruker Manual 2016 Bruker Manual 2005 Bruker 2019 Oct Service Training	Size Type Parent File fo File fo	Name L data smart Our NMR data Server disk
<	>	<
0 B of 22.3 MB in 0 of 28	1 hidden	0 B of 0 B in 0 of 2

Then you can download your data to your local disk.

3. Data Folder's location in each server:

<u>AV700, AV600, AV500 & AV300</u>: /home/opt/your group's NMR account name/ <u>NEO600 & NEO400</u>: /opt/ your group's NMR account name/ Usually it will have a few sub-folders and you can find your own data folder there.